

## Analyse mit Exomiser:

### 1. yml Datei erstellen mit Pfad zur ped, vcf, proband sample name und phänotypen

```
## Exomiser Analysis Template.
# These are all the possible options for running exomiser. Use this as a template for
# your own set-up.
---
analysis:
  ped: /PATH/TO/PED
  proband: PROAND_SAMPLE_NAME
  vcf: /PATH/TO/VCF
  # AUTOSOMAL_DOMINANT, AUTOSOMAL_RECESSIVE, X_RECESSIVE or UNDEFINED
  inheritanceModes: {
    AUTOSOMAL_DOMINANT: 0.1,
    AUTOSOMAL_RECESSIVE_HOM_ALT: 3.0,
    AUTOSOMAL_RECESSIVE_COMP_HET: 5.0,
    X_DOMINANT: 0.1,
    X_RECESSIVE_HOM_ALT: 3.0,
    X_RECESSIVE_COMP_HET: 5.0,
  }
  #FULL, SPARSE or PASS_ONLY
  analysisMode: PASS_ONLY
  #hpoIds: ['HP:0100280']
  hpoIds: []
  #Possible frequencySources:
  #Thousand Genomes project http://www.1000genomes.org/
  # THOUSAND_GENOMES,
  #ESP project http://evs.gs.washington.edu/EVS/
  # ESP_AFRICAN_AMERICAN, ESP_EUROPEAN_AMERICAN, ESP_ALL,
  #ExAC project http://exac.broadinstitute.org/about
  # EXAC_AFRICAN_INC_AFRICAN_AMERICAN, EXAC_AMERICAN,
  # EXAC_SOUTH_ASIAN, EXAC_EAST_ASIAN,
  # EXAC_FINNISH, EXAC_NON_FINNISH_EUROPEAN,
  # EXAC_OTHER
  frequencySources: [
    THOUSAND_GENOMES,
    TOPMED,
    UK10K,
    ESP_AFRICAN_AMERICAN, ESP_EUROPEAN_AMERICAN, ESP_ALL,
    EXAC_AFRICAN_INC_AFRICAN_AMERICAN, EXAC_AMERICAN,
    EXAC_SOUTH_ASIAN, EXAC_EAST_ASIAN,
    EXAC_FINNISH, EXAC_NON_FINNISH_EUROPEAN,
    EXAC_OTHER,
    GNOMAD_E_AFR,
    GNOMAD_E_AMR,
    # GNOMAD_E_ASJ,
    GNOMAD_E_EAS,
    GNOMAD_E_FIN,
    GNOMAD_E_NFE,
    GNOMAD_E_OTH,
    GNOMAD_E_SAS,
    GNOMAD_G_AFR,
    GNOMAD_G_AMR,
    # GNOMAD_G_ASJ,
    GNOMAD_G_EAS,
    GNOMAD_G_FIN,
    GNOMAD_G_NFE,
    GNOMAD_G_OTH,
    GNOMAD_G_SAS
  ]
  #Possible pathogenicitySources: POLYPHEN, MUTATION_TASTER, SIFT, CADD, REMM
  #*WARNING* if you enable CADD, ensure that you have downloaded and installed the CADD tabix files
  #and updated their location in the application.properties. Exomiser will not run without this.
  pathogenicitySources: [POLYPHEN, MUTATION_TASTER, SIFT]
  #this is the standard exomiser order.
  #all steps are optional
  steps: [
    #intervalFilter: {interval: 'chr10:123256200-123256300'},
    #genePanelFilter: {geneSymbols: []},
    #failedVariantFilter: {},
    qualityFilter: {minQuality: 20.0},
    variantEffectFilter: {remove: [UPSTREAM_GENE_VARIANT,
      INTERGENIC_VARIANT,
      NON_CODING_TRANSCRIPT_INTRON_VARIANT,
      SYNONYMOUS_VARIANT,
      DOWNSTREAM_GENE_VARIANT]},
    #knownVariantFilter: {}, #removes variants represented in the database
    frequencyFilter: {maxFrequency: 5.0},
    pathogenicityFilter: {keepNonPathogenic: true},
    #inheritanceFilter and omimPrioritiser should always run AFTER all other filters have completed
    #they will analyse genes according to the specified modeOfInheritance above- UNDEFINED will not be analysed.
    inheritanceFilter: {},
    #omimPrioritiser isn't mandatory.
    omimPrioritiser: {},
    #priorityScoreFilter: {minPriorityScore: 0.4},
    #Other prioritisers: Only combine omimPrioritiser with one of these.
    #Don't include any if you only want to filter the variants.
    hiPhivePrioritiser: {},
    # or run hiPhive in benchmarking mode:
```

```

        #hiPhivePrioritiser: {runParams: 'mouse'},
        #phivePrioritiser: {}
        #phenixPrioritiser: {}
        #exomeWalkerPrioritiser: {seedGeneIds: [11111, 22222, 33333]}
    ]
    outputOptions:
        outputPassVariantsOnly: false
        #numGenes options: 0 = all or specify a limit e.g. 500 for the first 500 results
        numGenes: 0
        #outputPrefix options: specify the path/filename without an extension and this will be added
        # according to the outputFormats option. If unspecified this will default to the following:
        # {exomiserDir}/results/input-vcf-name-exomiser-results.html
        # alternatively, specify a fully qualified path only. e.g. /users/jules/exomes/analysis
        outputPrefix: analysis/results
        #out-format options: HTML, TSV-GENE, TSV-VARIANT, VCF (default: HTML)
        outputFormats: [TSV-GENE, TSV-VARIANT, VCF, HTML]

```

## 2. Starten des Durchlaufs

Auf dem cluster mit:

```

module load IKMB Java/1.8.0 && sbatch --exclude=rzcl236 -p ikmb_a --qos ikmb_a -J exomiser --mem=16G -t 2:0:0 --
wrap="java -Xms8g -Xmx16g -jar /ifs/data/nfs_share/ikmb_repository/software/exomiser/10.0.0/exomiser-cli-10.0.0.jar
--analysis analysis/config_exomiser.yml"

```

## 3. Ausführen von run\_pipeline.py

Parameter:

```

exomiser_path = '/Users/broder/projects/IBD-DACH_FT_3/exomiser/analysis'
annovar_path = '/Users/broder/projects/IBD-DACH_FT_3/variants.merged.normalized.vcf.gz.annovar.hg19_multianno.txt'
output_path = '/Users/broder/projects/IBD-DACH_FT_3/result.tsv'

```

- Merged die Varianten von den verschiedenen Exomiser Inheritance Dateien zusammen

H
FILTER
AD_PASS;AR_PASS
AR_PASS
AR_PASS
AD_PASS
AR_PASS
AR_PASS
AR_PASS
AR_PASS
AR_PASS
AR_PASS
AR_PASS
AD_PASS;AR_PASS
AD_PASS;AR_PASS

- Wendet die verschiedenen Genlisten auf die Varianten an

GeneFilter
ClinVar;imm
pedibdb
pedibdb
ClinVar;DisGeNET
ClinVar;DisGeNET
ClinVar;imm
DisGeNET
DisGeNET
DisGeNET
DisGeNET
DisGeNET
DisGeNET
ClinVar;pedibdb
ClinVar;imm
ClinVar;imm
ClinVar;pedibdb;imm
ClinVar;pedibdb;imm
ClinVar;imm
ClinVar;imm

- Fügt die AnnoVar Annotation für jede Variante hinten an die Exomiser Outputdatei

#### 4. Kopieren der Ergebnisse in das Result Template

#### 5. Filtern

A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S		
keep	note	#CHROM	POS	REF	ALT	QUAL	FILTER	GENOTYP	COVERAGE	FUNCTION	HGV	EXOMISER	CADD	POLYPHE	MUTATIO	SIFT	REMM	DBSNP	IF	MA
1	https://www	chr7	87170711	C	T	1221.41	AR_PASS	1/1:0/1	46	missense_va	ABCB1:ENST1	ABCB1	15.54	0.063	.	0.077	.	.	.	.
1	chr19	4218002	C	T	13848.4	AR_PASS	1/1:0/1	416	missense_va	ANKRD24:EN	ANKRD24	34	0.996	.	.	.	.	.	rs771153272	0.0
1	https://www	chr7	99704283	G	A	2948.14	AR_PASS	1/1:0/1	99	missense_va	AP4M1:ENST	AP4M1	20.4	0.011	0.999955	0.363	.	.	rs574878636	0.3
1	https://acadi	chr7	117232642	A	G	6865.11	AR_PASS	1/1:0/1	208	missense_va	CFTR:ENSTOC	CFTR	21.9	0.976	0.997072	0.001	.	.	rs1800103	0.0
1	chr10	105193754	G	A	1501.04	AR_PASS	1/1:0/1	69	missense_va	PDCD11:ENS	PDCD11	28	0.694	0.947106	0.127	.	.	rs61760969	0.2	
1	chr10	97192237	G	A	2675.13	AR_PASS	1/1:0/1	108	missense_va	SORBS1:ENS	SORBS1	29.5	1	1	0	.	.	rs200179325	1.	
1	chr12	72300873	C	T	1363.64	AR_PASS	1/1:0/1	60	missense_va	TBC1D15:EN	TBC1D15	32	1	1	0.002	.	.	rs777708288	0.0	
1	chr19	3728612	G	A	9479.93	AR_PASS	1/1:0/1	308	missense_va	TJP3:ENSTOO	TJP3	25.1	0.983	1	0.258	.	.	rs367879011	0.0	
1	chr10	98155678	G	A	3298.97	AR_PASS	1/1:0/1	138	missense_va	TLL2:ENSTOO	TLL2	17.83	0.96	0.98092	0.035	.	.	rs41291628	2.	
1	chr7	100371474	GTT	ATTG	6216.48	AR_PASS	1/1:0/1	250	frameshift_e	ZAN:ENSTOO	ZAN	.	.	.	.	.	.	.	.	.

#### 6. Extrahieren der gefilterten Varianten zu einem vcf

```

Sub toVcf()

Dim sh As Worksheet
Dim rw As Range
Dim keep As Integer

outputFile = FreeFile
Open "/Users/broder/tmp/vcf_output.vcf" For Output As #outputFile
Print #outputFile, "#CHROM" & vbTab & "POS" & vbTab & "ID" & vbTab & "REF" & vbTab & "ALT" & vbTab & "QUAL" & vbTab & "FILTER" & vbTab & "INFO"

Set sh = ActiveSheet
start = 3
For Each rw In sh.Rows
    If sh.Cells(rw.Row, 1).Value = "1" Then
        chromTmp = sh.Cells(rw.Row, start).Value
        chrom = Replace(chromTmp, "chr", "")
        pos = sh.Cells(rw.Row, start + 1).Value
        id = sh.Cells(rw.Row, start + 16).Value
        ref = sh.Cells(rw.Row, start + 2).Value
        alt = sh.Cells(rw.Row, start + 3).Value
        qual = sh.Cells(rw.Row, start + 4).Value
        filterVal = sh.Cells(rw.Row, start + 5).Value
        info = "."
        Print #outputFile, chrom & vbTab & pos & vbTab & id & vbTab & ref & vbTab & alt & vbTab & qual & vbTab & filterVal & vbTab & info
    ElseIf sh.Cells(rw.Row, start).Value = "" Then
        Exit For
    End If
Next rw
Close #outputFile
End Sub

```

-> Einstellen in VarWatch